



# ***TA Instruments Q800 to DMA 850 Transition Guide***

---

## **Overview**

For users upgrading from the Q Series™ Q800 DMA to the new Discovery DMA 850 Dynamic Mechanical Analyzer, this guide outlines the key changes between the two instruments and provides details on the information needed for a smooth transition to the DMA 850.

## **Test Setup and Procedures**

Test procedures could be setup and accessed on the Q800 by selecting an operating mode such as "DMA Multi-strain" and then accessing the tests within that mode. In each operating mode, a number of tests were available - this included pre-programmed standard tests as well as the option to customize the test procedure through the Custom mode. The DMA 850 expands on this tradition of easy to use interfaces through the introduction of DMA Express and the DMA Unlimited procedures - two new and intuitive interfaces for user interaction.

DMA Express procedures contain all the pre-programmed tests previously available on the Q800 neatly packaged into single step procedures. This mode, ideal for performing simple tests that do not require customization, is designed specifically to get new users up and running on the DMA as fast as possible. Each test form has been optimized for maximum efficiency with all the key parameters needed to run the test presented in a logical and sequential format.

DMA Unlimited procedures expand on the fundamental test modes available in DMA Express by allowing access to additional conditioning steps and the ability customize test methods. In addition, it also introduces the new capability to sequence DMA procedures that can mix and match different operating modes within a single test. For example, tests containing both oscillation and transient (creep, stress relaxation) steps can now be programmed on the DMA 850 thanks to this advanced mode of operation, a capability that was not available with previous generations.

The fundamental test blocks that are needed to program the test are common to both procedures - all the tests are laid out in one of four operating modes. A summary of the available test modes in DMA Express and DMA Unlimited are given below:

- Oscillation
  - 1 Frequency Sweep
  - 2 Strain Sweep
  - 3 Stress Sweep
  - 4 Temperature Sweep

- 5 Temperature Ramp
- 6 Time Sweep
- 7 Temperature Sweep (Multifrequency)
- 8 Fatigue Test
- 9 Temperature Ramp (Multifrequency)
- 10 Temperature Ramp (Multistep)
- Strain Control
  - 1 Stress Relaxation
  - 2 Stress Relaxation TTS
  - 3 IsoStrain
- Stress Control
  - 1 Creep
  - 2 Creep Recovery
  - 3 Creep TTS
  - 4 IsoStress
- Rate Control
  - 1 Strain Ramp
  - 2 Stress Ramp

If you are upgrading to the Discovery DMA 850 from the Q Series, use the information presented in the table below to determine the DMA 850 mode that is equivalent to the tests used for the Q800.

## Standard DMA Tests–DMA Express, DMA Unlimited.

| Q800 Test                  |                            | Equivalent DMA 850 Test |                  |                           | Comments         |
|----------------------------|----------------------------|-------------------------|------------------|---------------------------|------------------|
| Mode                       | Test                       | Mode                    | Test             | Controlled Test Parameter |                  |
| DMA Multi-Strain           | Strain Sweep               | Oscillation             | Strain Sweep     | Amplitude or Strain       |                  |
| DMA Multi-Stress           | Stress Sweep               | Oscillation             | Stress Sweep     | Force or Stress           |                  |
| DMA Multi-Frequency Strain | Isothermal Temp/Freq Sweep | Oscillation             | Frequency Sweep  | Amplitude or Strain       |                  |
| DMA Multi-Frequency Strain | Temp Ramp/Freq Sweep       | Oscillation             | Temperature Ramp | Amplitude or Strain       | Single frequency |

| Q800 Test                  |                            | Equivalent DMA 850 Test |                                       |                                     | Comments   |
|----------------------------|----------------------------|-------------------------|---------------------------------------|-------------------------------------|--|
| Mode                       | Test                       | Mode                    | Test                                  | Controlled Test Parameter           |  |
| DMA Multi-Frequency Strain | Temp Ramp/<br>Freq Sweep   | Oscillation             | Temperature Ramp<br>(Multifrequency)  | Amplitude or Strain                 | Multiple frequencies   |
| DMA Multi-Frequency Strain | Temp Ramp/<br>Freq Sweep   | Oscillation             | Temperature Sweep                     | Amplitude or Strain                 | Single frequency   |
| DMA Multi-Frequency Strain | Temp Ramp/<br>Freq Sweep   | Oscillation             | Temperature Sweep<br>(Multifrequency) | Amplitude or Strain                 | Multiple frequencies   |
| DMA Multi-Frequency Stress | Isothermal Temp/Freq Sweep | Oscillation             | Frequency Sweep                       | Force or Stress                     |  |
| DMA Multi-Frequency Stress | Temp Ramp/<br>Freq Sweep   | Oscillation             | Temperature Ramp                      | Force or Stress                     | Single frequency   |
| DMA Multi-Frequency Stress | Temp Ramp/<br>Freq Sweep   | Oscillation             | Temperature Ramp<br>(Multifrequency)  | Force or Stress                     | Multiple frequencies   |
| DMA Multi-Frequency Stress | Temp Step/<br>Freq Sweep   | Oscillation             | Temperature Sweep                     | Force or Stress                     | Single frequency   |
| DMA Multi-Frequency Stress | Temp Step/<br>Freq Sweep   | Oscillation             | Temperature Sweep<br>(Multifrequency) | Force or Stress                     | Multiple frequencies   |
| N/A                        |                            | Oscillation             | Time Sweep                            | Amplitude, Strain, Force, or Stress | New test, Oscillation test at isothermal temperature, constant test conditions   |
|                            |                            | Oscillation             | Fatigue                               | Amplitude, Strain, Force, or Stress | New test, Oscillation test at isothermal temperature, constant test conditions   |
|                            |                            | Oscillation             | Temperature Ramp<br>(MultiStep)       | Amplitude, Strain, Force, or Stress | New test, DMA Express only: Allows the sequencing of multiple temperature ramps in the same test, such as heat-cool-heat |
| DMA Creep                  | Creep                      | Stress Control          | Creep                                 | Force or Stress                     | Only creep step, no recovery   |

| Q800 Test             |                                | Equivalent DMA 850 Test |                       |                           | Comments                                     |
|-----------------------|--------------------------------|-------------------------|-----------------------|---------------------------|--|
| Mode                  | Test                           | Mode                    | Test                  | Controlled Test Parameter |  |
| DMA Creep             | Creep                          | Stress Control          | Creep Recovery        | Force or Stress           | Combines creep and recovery in a single step |
| DMA Creep             | Creep TTS                      | Stress Control          | Creep TTS             | Force or Stress           |  |
| DMA Controlled Force  | Temp Ramp/<br>Controlled Force | Stress Control          | Isostress             | Force or Stress           |  |
| DMA Stress Relaxation | Stress Relaxation              | Strain Control          | Stress Relaxation     | Displacement or Strain    |  |
| DMA Stress Relaxation | Stress Relaxation TTS          | Strain Control          | Stress Relaxation TTS | Displacement or Strain    |  |
| DMA Iso-Strain        | Isostrain                      | Strain Control          | IsoStrain             | Displacement or Strain    |  |
| DMA Controlled Force  | Stress/Strain                  | Rate Control            | Stress Ramp           | Force or Stress           |  |
| DMA Strain Rate       | Strain Ramp                    | Rate Control            | Strain Ramp           | Displacement or Strain    |  |

## Customized DMA Tests & Sample Conditioning – DMA Unlimited

In addition to the standard DMA tests, DMA Unlimited procedures also have sample conditioning blocks that be sequenced together to generate customized multi-step procedures similar to the Custom tests available within each mode on the Q800. All basic DMA tests – Oscillation, Strain Control, Stress Control, Rate Control – are available in this mode as well.

Custom tests on the Q800 were limited to the selected DMA mode only - however, this is not a limitation with the DMA 850 since the DMA Unlimited procedures support the sequencing of tests from multiple operating modes within a single procedure. The individual segments available through the Custom tests on the Q800 can be found within the Conditioning mode on the DMA 850.

The information presented in the table below summarizes the DMA 850 modes that are equivalent to segments available on the Q800.

| Q800 Custom Segments  | Equivalent DMA 850 Test |             |   |   |
|-----------------------|-------------------------|-------------|---|---|
| Segment               | Mode                    | Test        | Option/<br>Controlled Test<br>Parameter | Comments  |
| Jump                  | Conditioning            | Temperature | Jump                                    |   |
| Equilibrate           | Conditioning            | Temperature | Equilibrate                             |   |
| Initial temperature   | Conditioning            | Temperature | Initial                                 |   |
| Ramp (temperature)    | N/A                     |             |   | Ability to ramp temperature is available within individual tests in each control mode                         |
| Isothermal            | N/A                     |             |   | Ability to perform isothermal tests is available within individual tests in each control mode                 |
| Increment temperature | Conditioning            | Repeat      |   |   |
| Repeat                | Conditioning            | Repeat      |   |   |
| Repeat until          | N/A                     |             |   | Repeat blocks can be terminated through Limit Checking available within individual tests in each control mode |
| Motor drive           | Conditioning            | Other       | Motor                                   |   |
| Step (temperature)    | N/A                     |             |   | Ability to perform temperature steps (sweeps) is available within individual tests in each control mode       |
| Abort next segment    | N/A                     |             |   | Steps can be terminated through Limit Checking available within individual tests in each control mode         |
| Sample interval       | N/A                     |             |   | Data sampling intervals are directly set within individual tests in each control mode                         |
| External event        | Conditioning            | Other       | External Event                          |   |
| Data Storage          | Conditioning            | Data        | Data Storage                            |   |
| Fill cooler           | Conditioning            | GCA         |   |   |
| Heater PID            | Conditioning            | Other       | Override Heater PID                     |   |

| Q800 Custom Segments | Equivalent DMA 850 Test |                 |   |          |
|----------------------|-------------------------|-----------------|---|----------|
| Segment              | Mode                    | Test            | Option/<br>Controlled Test<br>Parameter | Comments |
| Frequency sweep      | Oscillation             | Frequency sweep | Amplitude, Strain, Force, or Stress     |          |
| Strain sweep         | Oscillation             | Strain Sweep    | Amplitude or Strain                     |          |
| Displace/Recover     | Stress Control          | Creep Recovery  | Force or Stress                         |          |
| Force                | Conditioning            | Stress          | Force or Stress                         |          |
| Ramp force           | Rate Control            | Stress Ramp     | Force or Stress                         |          |
| Increment force      | Conditioning            | Stress          | Force or Stress                         |          |
| Isostrain            | Strain Control          | IsoStrain       | Displacement or Strain                  |          |
| Set strain           | Conditioning            | Strain          | Displacement or Strain                  |          |
| Ramp strain          | Rate Control            | Strain Ramp     | Displacement or Strain                  |          |
| Ramp displacement    | Rate Control            | Strain Ramp     | Displacement or Strain                  |          |
| Stress               | Conditioning            | Stress          | Force or Stress                         |          |
| Ramp Stress          | Rate Control            | Stress Ramp     | Force or Stress                         |          |
| Humidity             | Conditioning            | Humidity        | Set RH                                  |          |
| Increment humidity   | Conditioning            | Humidity        | Increment RH                            |          |
| Step humidity        |                         |                 |   |          |

## Importing Q800 Data and Signal Mapping

The capability to natively import Q800 data is available within TRIOS versions 4.3 or later. When a Q800 file is opened in TRIOS, a large number of signals are presented with the same variable name as the Q800.

However, a limited number of signals are mapped to the newer DMA 850 variable names in the interest of consistency between the two instruments. This provides the ability to directly overlay Q800 data with DMA 850 data within TRIOS. Note that the DMA 850 data can also be exported and opened within Universal Analysis for backward compatibility.

The information presented in the table below summarizes the Q800 signals and the corresponding name used for imported data in TRIOS.

## Common Signals

| Q800 Signal | Imported Name in TRIOS | Units    | Comments   |
|-------------|------------------------|----------|--|
| Time        | Time/Step time         | s or min | Time and step time are the same in imported Q800 data; the DMA 850 reports individual step times in the DMA Unlimited mode |
| Temperature | Temperature            | °C       |  |

Signals not listed above are imported and mapped to the same name as the Q800 - examples include Storage Modulus, Loss Modulus,  $\tan(\delta)$ , etc.

## Oscillation Signals

| Q800 Signal         | Imported Name in TRIOS    | Units   | Comments  |
|---------------------|---------------------------|---------|---|
| Drive force         | Oscillation Force (drive) | N       | Force amplitude applied to the DMA motor to produce the commanded oscillation displacement/ force   |
| Force/Dynamic Force | Oscillation Force         | N       | Force exerted by the clamp on the sample when oscillating at the commanded oscillation displacement/ force  |
| Static Force        | Axial Force               | N       | Absolute force applied by the DMA motor on the sample to maintain contact with the sample (compression) or hold the sample in tension (film)  |
| Displacement        | $\Delta L$                | mm      | Refers to the change in the drive position of the DMA shaft during the test   |
| Length              | Position                  | mm      | Refers to the position signal visible on the DMA 850 user interface- this signal corresponds to the sample length in tension and the sample thickness in compression testing (see note below) |
| Position            | Drive-position            | mm      | Refers to the absolute position of the DMA drive shaft (see note below)   |
| Amplitude           | Oscillation Displacement  | $\mu m$ |   |
| Stress              | Oscillation Stress        | N       |   |
| Strain              | Oscillation Strain        | %       |   |
| N/A                 | Pretension ration         | %       | Refers to the ratio of the axial force to the oscillation force   |

Signals not listed above are imported and mapped to the same name as the Q800 - examples include Storage Modulus, Loss Modulus,  $\tan(\delta)$ , etc.

## Non-Oscillation Signals (All transient tests such as Creep, Stress Relaxation, IsoStrain, etc.)

| Q800 Signal        | Imported Name in TRIOS | Units         | Comments  |
|--------------------|------------------------|---------------|---|
| Static Force       | Total Force            | N             | Total force applied on the sample   |
| N/A                | Force                  | N             | Refers to the force applied on the sample during the current step, after subtracting the preload force (static force)   |
| Displacement       | Displacement           | $\mu\text{m}$ | (see note below)  |
| Position           | Drive position         |               | Refers to the absolute position of the DMA drive shaft (see note below)   |
| Length             | Position/<br>Length    | mm            | Refers to the position signal visible on the DMA 850 user interface—a separate length signal is also available in film, compression clamp data (see note below) |
| Strain             | Strain                 | %             |   |
| Stress             | Stress                 | MPa           | Refers to the additional stress on the sample after subtracting the preload force— Creep, Stress Relaxation, IsoStrain, IsoStress tests only                    |
| Stress             | Total Stress           | MPa           | Refers to the total stress on the sample, including the preload force - Stress Ramp, Strain Ramp only   |
| Strain Recovery    | Recovered Strain       | %             | Only available in creep tests   |
| Relaxation Modulus | Modulus                | MPa           | Available in all transient tests except creep   |

### Notes (all variable names refer to the Imported Name in TRIOS):

- 1 The Q800 reported a zero value for the Position signal in all clamps except film/fiber tension and compression. In such cases (cantilever, 3-point bending, shear sandwich), it is recommended to plot the  $\Delta L$  signal instead.
- 2 The Q800's Drive position reference was 0 mm at the top of travel, while the DMA 850's 0 mm reference is at the bottom of the travel. While overlaying data from the Q800, the trend for this signal will likely be opposite to that of the DMA 850. It is recommended to plot the Position or the Length signal in such cases.
- 3 For compression clamp data, the Displacement signal is negative on the Q800 while the DMA 850 uses positive notation for displacement in all tests. To overlay data, it is recommended to plot the absolute values of the displacement signal in TRIOS.