



ARES-G2 / RSA-G2 LN2 Kit Installation Instructions

Parts Supplied

Table 1 Parts supplied

Description	Quantity
LN2 Installation Instructions	1
LN2 Declaration of Conformity	1
LN2 cooler assembly	1
LN2 transfer line	1
LN2 flex line	1
Tube	1
Brass adapter 1/4-20 NPT	1
Coupling straight 1/4 NPT brass	2
Conn m thermocouple 1/4 NPT 3/8 tb SST	2
Relief valve 13 psi	1

Tools Required

- Tool kit included with ARES-G2 or RSA-G2 instrument
- Adjustable wrench for properly tightening fittings

Installation Procedure

This procedure explains how to install the LN2 Accessory Kit on the TA Instruments ARES-G2 or RSA-G2. The following procedure assumes an FCO mounted on the instrument.

Prepare the Instrument

- 1** Turn off the instrument and unplug the power cord.
- 2** Turn off all gas supplies at their source.

Prepare the LN2 Cooler

- 1 Remove the LN2 Cooler front panel by removing the four screws holding it in place.
- 2 From the LN2 ship kit, remove the LN2 flex line and relief valve assembly. Refer to [Figure 3](#) for part identification, if necessary.



LN2 flex line

Figure 1 Flex line shown.

- 3 Using the adjustable wrench, remove the fitting and plunger from the center NPT connection in the LN2 Cooler. See the figure below.



Fitting and plunger
to remove

Figure 2 Fitting and plunger in LN2 Cooler.

- 4 Using the parts obtained in [step 2](#) above, install the parts in the configuration shown below.




Figure 3 Tube, adapter, couplings, thermocouples, and relief valve installed in proper locations.

- 5 Replace the LN2 Cooler front panel using the four previously-removed screws and proceed to the next section.

Connect the LN2 Cooler to the Instrument

- 1 Connect one end of the LN2 flex line to the LN2 Cooler supply inlet filter connection (location shown below in [Figure 4](#)) and the other end to a liquid nitrogen supply source (bulk tank).

CAUTION: The bulk tank must be a low pressure (22 psi / 1.5 bars) source, as a high pressure tank will damage the LN2 controller. Verify that the tank is a low pressure source before completing step 1.



MISE EN GARDE: Le réservoir en vrac doit être une source de basse pression (22 psi / 1,5 bar), vu qu'un réservoir haute pression endommagerait le régulateur LN2. Assurez-vous que le réservoir est une source de basse pression avant de réaliser l'étape 1.




Figure 4 LN2 flex line and Cooler connection.

- 2 Using the adjustable wrench, remove the pipe cap from the FCO LN2 connection (circled below in [Figure 5](#)).




Figure 5 LN2 connection and pipe cap location (circled) on back of right-hand mounted FCO.

- 3 Connect the **elbowed end** of the LN2 transfer line (the hose with the larger diameter) to the FCO LN2 supply connection. Refer to [Figure 8](#) for left-hand and right-hand mounted FCO configurations.




Figure 6 Transfer line connected to FCO LN2 supply connection (right-hand mounted FCO shown).

- 4 Connect the opposite end of the transfer line to the LN2 Cooler, as shown below.




Figure 7 LN2 transfer line connection on LN2 Cooler.

The configuration should resemble the figures below.




Figure 8 LN2 Cooler configuration with left-hand and right-hand mounted FCO.

- 5 Plug the LN2 Cooler control cable into the LN2 connector on the FCO Controller.




Figure 9 LN2 control cable connection location on FCO Controller.

Restore the Instrument to Operating Condition

- 1 Turn on gas supplies and verify correct air pressures.
- 2 Plug in and power on the instrument.



NOTE: When the LN2 installation is complete, be sure to open the liquid nitrogen tank's liquid output valve before operating the LN2.
