

Application & Instrument Training TRIBOMETER

@ ITESM - CEM
Laboratorio de Tribología



Gustavo Ruiz
&
Osamu Katagiri

Introduction

Discovery Hybrid Rheometer (DHR) from TA Instruments.

- Optical Encoder Dual Reader.
 - Resolution better than 2 nanoradians.
- Advanced drag cup motor.
 - Keeps inertia, temperature and friction an absolute minimum.
- 2nd generation magnetic bearing.
 - Less than 70% of the friction of air-bearing.
 - It can measure down to 0.5 nNm of torque.
- Force rebalance transducer (FRT).
 - High accuracy normal force measurements up to 50N.
- True Position Sensor (TPS).
 - True gap accuracy via a high resolution linear position sensor.



Figure 1. Discovery hybrid Rheometer. TA Instruments.



Purposes

- To learn how to turn on the instrument in the correct order and the conditions needed for using it.
- To understand the software and how to design Rheology experiments.
- To learn how to calibrate the instrument before each experiment.
- To perform rheology experiments for practicing.



Materials




Certified Viscosity Standard. S600.



Discovery hybrid Rheometer.
TA Instruments.



Peltier plate for rheology.



Geometry: Parallel plate 60 mm.



Ethanol.



Procedure

➤ Initial setup



remove
protective covers



open air
valves
AKA. ball
valve / T




check air pressure, which
shall be around 30 (± 2) psi

SW UP
is ON

DOWN
is OFF



turn on the TA-
supplied Air
Cooled
Circulator (below
the working table)



plug the
instruments
to the orange
(isolated ground)
outlet



Procedure

➤ Initial setup (continued)



remove the
bearing clamp



turn on the
Electronics
Control
Module



"I"
is ON

"O"
is OFF

➤ Replace the Lower Stage (if required)

Press the
Release button



disconnect the
power cable



disconnect
the
fluid hoses



Procedure

➤ Replace the Lower Stage (continued)

Press the
Release button
again



remove the
attachment



Press the
Release button
again



align and fit
the attachment



connect the
power cable
& fluid hoses



Procedure

➤ Install a Geometry



clean the
geometry
(with alcohol)




The
geometry is
gently placed
(clockwise
direction
until the first
stop)



Procedure


➤ Prepare the sample



the sample is placed on the Peltier plate



the geometry is lowered using the arrows (until it is touching the sample)



the sample that comes out of the geometry, it is cleaned



It is checked that all the area of the geometry is covered by the sample




Procedure


➤ Run the experiment



Open the TA
software &
connect to the
instrument




input the experiment
data




select the geometry



calibrate the
instrument's
inertia




calibrate the
instrument's
friction




Procedure


➤ Run the experiment (continued)




Set the plate's temperature



perform the rotational mapping




Get the Zero Gap



run the experiment

define the experimental procedure




Results

Rheology Oil STD 25°C
Viscosity Measurement
AOKT
[Experiment 2]

Discovery HR-3, 12/03/2019 11:23:01 a. m.

60mm 0.9969" cone plate, Peltier plate Steel

Rheology Oil STD 25°C



$$\ln[1] = \text{error} = \frac{\text{Abs}[1032.71 - 1037]}{1037} \times 100$$

Out[1] = 0.413693

CANNON
INSTRUMENT COMPANY

2139 High Tech Road
State College, PA 16803
814-353-8000 • 800-476-6232 • Fax 814-353-9007
cannon@cannoninstrument.com
www.cannoninstrument.com

CERTIFICATE OF ANALYSIS

CANNON® CERTIFIED VISCOSITY REFERENCE STANDARD					
Viscosity Standard: S600			Lot Number: 18101		
Certification/Issue Date: 02/27/2018			Expiry Date: 02/28/2020		
Temperature ° C	° F	Kinematic Viscosity mm ² /s (cSt)	Dynamic Viscosity mPa.s (cP)	Density g/cm ³ (g/mL)	Saybolt Viscosity seconds
20.00	68.00	1718	1454	0.8466	145 SFS
25.00	77.00	1229	1037	0.8436	
37.78	100.00	575.5	481.1	0.8360	
40.00	104.00	509.8	425.5	0.8347	
50.00	122.00	308.0	255.3	0.8287	
80.00	176.00	91.08	73.87	0.8111	236 SUS
98.89	210.00	50.50	40.40	0.8000	
100.00	212.00	48.87	39.06	0.7993	

Tested and certified in the U.S.A.

